

TECHSPEC[®] 125mm, Center Drive, Thru-Hole, Metric Stage



TECHSPEC[®] Crossed-Roller Bearing Linear Translation Stages (Thru-Hole Model)



Stock #66-468 3 In Stock

-

1

+

S\$994⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-2	S\$994.00 each
Qty 3+	S\$867.30 each
Need More?	Request Quote

Product Downloads

General

Type:

Metric

Physical & Mechanical Properties

Type of Movement:

Linear (X)

Guide System:



