

[See all 6 Products in Family](#)

3" Length Bag, Cotton-Tipped Swabs



3" Length Bag, Cotton-Tipped Swabs (1 Bag)

Stock **#56-926** **1 In Stock**

-

1

+

S\$8⁴⁰

ADD TO CART

Volume Pricing	
Qty 1+	S\$8.40 each
Need More?	Request Quote

Product Downloads

General	
1 Bag	Type:
100	Number of Swabs:
Physical & Mechanical Properties	
3.00	Length (inches):
Regulatory Compliance	

Product Details

- Lint-Free
- 3" or 6" Wooden Stems
- Bags, Boxes, and Cases Available

The Cotton-Tipped Swab Applicators are ideal for optical cleanup applications. They are offered with either 3" or 6" wooden stems. These non-sterile swabs are great for providing versatility in a multitude of applications. Cotton-Tipped Swab Applicators are lint-free to be better functional in cleaning applications. Each bag contains 100 swabs, each box contains 1,000 swabs, and each case contains 10,000 swabs.