

[See all 6 Products in Family](#)

## 3" Length Case (10 Boxes), Cotton-Tipped Swabs



3" Length Case (10 Boxes), Cotton-Tipped Swabs

Stock #56-928 **12 In Stock**

⊖ 1 ⊕ **\$221.20**

**ADD TO CART**

| Volume Pricing |                               |
|----------------|-------------------------------|
| Qty 1+         | <b>\$221.20</b> each          |
| Need More?     | <a href="#">Request Quote</a> |

### Product Downloads

#### General

Case (10 Boxes, 10 Bags per Box) **Type:**  
10,000 **Number of Swabs:**

#### Physical & Mechanical Properties

3.00 **Length (inches):**

#### Regulatory Compliance

## Product Details

- Lint-Free
- 3" or 6" Wooden Stems
- Bags, Boxes, and Cases Available

The Cotton-Tipped Swab Applicators are ideal for optical cleanup applications. They are offered with either 3" or 6" wooden stems. These non-sterile swabs are great for providing versatility in a multitude of applications. Cotton-Tipped Swab Applicators are lint-free to be better functional in cleaning applications. Each bag contains 100 swabs, each box contains 1,000 swabs, and each case contains 10,000 swabs.

---

;