


TECHSPEC® 70mm, Center Drive, Thru-Hole, English Extended Micrometer



Stock #66-440 **2 In Stock**

⊖ 1 ⊕ **\$896^{.00}**

ADD TO CART

Volume Pricing	
Qty 1-2	\$896.00 each
Qty 3+	\$781.48 each
Need More?	Request Quote

Product Downloads

General

English **Type:**

Physical & Mechanical Properties

Linear (X) **Type of Movement:**

Crossed Roller Bearing **Guide System:**

Center Drive **Drive Orientation:**

