


Replacement Tips for Resin 1.2mm Tip



Resin 1.2mm Tip Width, Non-Marring Tweezers, #36-449



Stock **#36-450** 3-5 DAYS

- 1 +

S\$21.75

ADD TO CART

Qty 1+
S\$21.75

Volume Pricing
Request Quote

Get Product Downloads

Length (mm): 48.00
Type: Resin Tip Replacement
Tip Width (mm): 1.2

Regulatory Compliance

RoHS:
[Compliant](#)

Certificate of Conformance:
[View](#)
