

TECHSPEC[®] Aperture Kit for Achromat Prototyping Kit



Aperture Kit for Achromat Prototyping Kit

Stock **#54-720** **1 In Stock**

-

1

+

S\$235²⁰

ADD TO CART

Volume Pricing	
Qty 1-5	S\$235.20 each
Qty 6-10	S\$218.82 each
Qty 11-25	S\$208.74 each
Need More?	Request Quote

Product Downloads

General	
Type:	Prototyping Kit
Function:	Kit

Regulatory Compliance	
RoHS 2015:	

