

TECHSPEC® DAPI (High Brightness) Filter Cube Set - Nikon



Olympus Filter Cube #86-833 with Pre-Mounted Fluorescence Filters



Stock #86-370-NKN [CONTACT US](#)

⊖ 1 ⊕ **S\$1,682.00**

ADD TO CART

Qty 1-9
S\$1,682.00

Qty 10+
S\$1,529.75

Volume Pricing
[Request Quote](#)

Get Product Downloads

Substrate:
[Fused Silica](#) (Coming 7980)

Coating:
Hard Coated

Type:
Fluorescence Filter Kit

Excitation Wavelength (nm):
352 - 402

Compatible Fluorophore:
DAPI (High Brightness)

Emission Wavelength (nm):
417 - 477

Dichroic Cut-On Wavelength (nm):
409.00

Excitation Filter:
[#84-093](#)

Emission Filter:
[#84-095](#)

Dichroic Filter:
[#86-330](#)

Wavelength Range (nm):
352 - 477

Compatible Microscope Brand:
Nikon

Regulatory Compliance

RoHS:
Not Compliant

Certificate of Conformance:
[View](#)
