

Diffuse reflectance phantom Brain 25 x 25 x 10 mm



Stock **#71-277** 1 In Stock

-

1

+

S\$945^{.00}

ADD TO CART

Volume Pricing	
Qty 1-4	S\$945.00 each
Qty 5-9	S\$850.50 each
Need More?	Request Quote

Product Downloads

SPECIFICATIONS

General

Type:

Brain

Physical & Mechanical Properties

Dimensions (mm):

25 x25 x 10

Material Properties

12Reduced Scattering Coefficient, μ_s' :

0.2Absorption Coefficient, μ_a :

Regulatory Compliance

[View](#)Certificate of Conformance:

PRODUCT DETAILS

- Mmic Muscle, Brain, or Adipose Tissue
- Available in 25 x25 mm and 50 x50 mm Sizes
- Ideal for a Range of Life Science Applications

Diffuse Reflectance Phantoms are solid samples that mimic the optical properties of biological tissues for testing, calibration, and characterization of biophotonics devices. These homogeneous phantoms, with well characterized optical properties at 780nm, simulate the optical properties and mechanical characteristics of biological tissues such as muscle, brain, and adipose tissue. These samples come in two sizes, 25 x25mm or 50 x 50mm, to accommodate multiple experimental setups. Diffuse Reflectance Phantoms are ideal for spatial frequency domain imaging, diffuse reflectance spectroscopy, diffuse Raman spectroscopy, and biophotonics/biomedical imaging.