

[See all 5 Products in Family](#)

Olympus MPLN 5X Objective

See More by [Olympus](#)



Stock #87-014 NEW **1 In Stock**

⊖ 1 ⊕ **SS812⁰⁰**

ADD TO CART

Volume Pricing

Qty 1+	SS812.00 each
Need More?	Request Quote

Product Downloads

General

Model Number:
MPLN5X

Compatible Tube Lens Focal Length (mm):
Focal Length: 180mm

Type:
Microscope Objective

Style:
Infinity Corrected

Manufacturer:
Olympus

Physical & Mechanical Properties

4.4 Field of View (mm):

25.00 Length excluding Threads (mm):

24 Maximum Diameter (mm):

64 Weight (g):

Optical Properties

N/A Compatible Cover Glass Thickness (mm):

36.00 Focal Length FL (mm):

5X Magnification:

0.10 Numerical Aperture NA:

3.35 Resolving Power (μm):

27.50 Depth of Field (μm):

20 Working Distance (mm):

400 - 700 Wavelength Range (nm):

22 Field Number (mm):

45 Parfocal Length (mm):

N/A Immersion Liquid:

7.20 Entrance Pupil Diameter (mm):

Threading & Mounting

RMS / 20.32mm x 36 TPI Mounting Threads:

Regulatory Compliance

[View](#) Certificate of Conformance:

Product Details

- Ideal for Transmitted or Reflected Brightfield Microscopy
- Produce Flat Images Up to Field Number (F.N.) 22
- Excellent Price-to-Performance Ratio
- Additional [Olympus Infinity Corrected Objectives](#) Available

Olympus MPlan Achromatic Objectives are designed for brightfield observations, providing excellent contrast and optimum flatness throughout the field of view. Engineered for inspection of silicon wafers, metals, and other industrial microscopy applications, these objectives may also be used in Raman microscopes and other reflected light techniques without coverslip correction. Olympus MPlan Achromatic Objectives are available in 5X-100X magnifications and provide excellent image flatness up to F.N. 22. [Olympus Plan Achromatic Objectives](#) are also available for observing biological specimens.