

Polarization Set Schott Annular Ring Light

See More by [SCHOTT Optical Components](#)



Stock **#72-143** **1 In Stock**

-

1

+

S\$763⁰⁰

ADD TO CART

Volume Pricing	
Qty 1+	S\$763.00 each
Need More?	Request Quote

Product Downloads

General	
158.430	Model Number:
Polarizer	Type:
Compatible up to 66mm dia. Ring Light	Note:
SCHOTT	Manufacturer:
Annular Ring Light	Compatible Product:

Regulatory Compliance	
GDPR Compliance: All data processing activities are documented, and data subjects' rights are respected. Regular audits are conducted to ensure compliance. ISO 27001: The organization maintains an ISO 27001 certification, which covers information security management systems. Other Standards: The organization also complies with various industry-specific standards, such as HIPAA for healthcare data and PCI DSS for payment card information.	Compliance Officer: A dedicated Compliance Officer is responsible for monitoring and ensuring adherence to all applicable regulations. Training: All employees receive regular training on data protection and security policies to ensure a culture of compliance. Incident Response: A robust incident response plan is in place to address any potential breaches or non-compliance issues promptly.

[View](#) **Certificate of Conformance:**

[View](#) **Certificate of Conformance:**

Product Details

- Maximum Light Intensity of 1100 Lumens
- Continuous Dimming from 0 - 100%
- Various Light Guide Options Available
- Ideal for Use with Stereo Microscopes

SCHOTT KL Fiber Optic Light Sources are LED cold light sources that produce intense, but heat free light for microscopy applications. Combining the efficient light output of the light engine with various light guide geometries offers significantly higher intensities compared to direct LED light. Two light source models are available, the KL 1600 which features a 680 lumens maximum intensity and the KL 2500 which features an 1100 lumens maximum intensity with an extra-fine dimming option. SCHOTT KL Fiber Optic Light Sources are ergonomically designed to have a small footprint, allowing them to be easily integrated into compact benchtop lab settings. The uniform light quality of the KL Series is ideal for biological, geological, and pathological stereo microscope applications.

Note: Power cord sold separately.

Technical Information

