

[See all 22 Products in Family](#)

Tube Lens



Stock #18-505 **12 In Stock**

− 1 + S\$728⁰⁰

ADD TO CART

Volume Pricing

Qty 1-4	S\$728.00 each
Qty 5+	S\$691.60 each
Need More?	Request Quote

Product Downloads

General

Note:
One required between camera tube and tube lens spacer/filter cube

Physical & Mechanical Properties

Length (mm):
16.15

Maximum Diameter (mm):

Optical Properties

120.00 Focal Length FL (mm):

1X Magnification:

Threading & Mounting

M29 x 0.75 & Ø30mm Dovetail Mounting Threads:

Regulatory Compliance

[View](#) Certificate of Conformance:

Product Details


- Up to 43% Smaller than Industry Standard Systems
- Optimized Performance up to 1.1" Sensor
- Ideal for Life Science and Industrial Applications
- Winner of the [2023 Bronze VSD Innovators Award](#)

Edmund Optics® 120i Plan APO Infinity Corrected Objectives are infinite conjugate microscope objectives that offer exceptional image flatness and chromatic correction over 400-700nm. These objectives were designed for a 1.1" sensor to offer a large field of view and reliable image quality for high-throughput imaging applications. Ideal for life science and machine vision applications such as microscopy, flow cytometry, pharmaceutical inspection, and assembly line or fault inspection.

These objectives won the [2023 Bronze VSD Innovators Award](#).

Technical Information

Compact Tube Lens System



Standard Tube Lens System

